

# D4000™

## Arthroscopic Resection System

# Features and Benefits

*ConMed's next generation of powered instrument consoles for arthroscopic resection, the **NEW D4000** has intuitive features that deliver new levels of simplicity, versatility and performance.*

### Easy-to-Use Interface

Designed to provide intuitive functionality, the D4000 dual touchscreen interface allows users to quickly and simply adjust a wide range of settings. Now adjusting the oscillation settings happens with the touch of a finger – eliminating inconvenient and time-consuming DIP switches.

### Automatic Fluid Adjustment

For a seamless clinical experience, the D4000 is directly integrated with ConMed's 24K® Pump. This means that these two units will effortlessly work together to automatically detect and adjust fluid levels – helping to maintain the surgical site and aid clinical visualization.



### Plug & Play Performance

Featuring two plug and play handpiece ports, the D4000 provides continuous, uninterrupted power. This console provides peak performance for both handpieces – even during simultaneous use at maximum RPMs.

### OR Versatility

A versatile console, the D4000 supports 31 handpieces, 5 languages, features built-in wireless connectivity to support multiple wireless foot controllers and allows for customizable user-profiles that let users save their preferred settings.

## Ordering Information

Description	Catalog Number
D4000 Arthro Power Console	D4000
24K Interface Cable	D4024

To schedule an evaluation, contact your local ConMed Sales Representative or call **1-800-237-0169**. To learn more about ConMed's resection products including blades, burs, handpieces and pumps, visit [ConMed.com](http://ConMed.com).

Local 727-392-6464  
Toll Free 1-800-237-0169

ConMed.com  
customer\_service@conmed.com